

Ponstan

Scan	Sensor Number
P2	6130007

Model: Ponstan, version 46cc79f86c2c4a12b8ce96ea9c066f80

Tablet = Bimastan

Ponstan

This model demonstrates a classification model for pharmaceutical tablets, using near-infrared (NIR) spectroscopy data in the 900nm - 1700nm wavelength range. The following tablets were used to build the model. Both are formulations of mefenamic acid. Ponstan is one of the most commonly counterfeited brands of pharmaceutical products in Indonesia, and the Bimastan tablets used to build this model were re-packaged with Pfizer's Ponstan brand.

Active Ingredient	Brand	Marks	Tablet Manufacturer
Mefenamic Acid 250mg	Ponstan	P-D	Pfizer Indonesia, Jakarta, Indonesia
Mefenamic Acid 250mg	Bimastan	BMF	Bima Mitra Farma, Tangerang , Banten, Indonesia

Disclaimer

The models used to create this prediction have been built by Sagitto Ltd using its best endeavours. As new data becomes available, new model versions may be created to improve model accuracy, and therefore results with future models may differ from those made with the current models. While Sagitto Ltd and the model owner have used their best endeavours to provide accurate predictions, neither Sagitto Ltd nor the model owner provide any guarantee of their accuracy. Sagitto Ltd and the model owner accept no liability for decisions made as a consequence of using the predictions from these models. The authenticity of this report can be verified by scanning the QR code.



Print date/time: 17/08/2016, 3:30:52 PM